

Biological Spill Kit

Every laboratory using Biological Agent(s) must have easy access to a Biological Spill kit. Place the Biological Spill kit(s) near the Biological Safety Cabinet (BSC) so that they are easily accessible in the event of a spill. Spill kits must be checked periodically, and restored after each use.

Items needed as a minimum are:

Gloves (2 pairs)

Safety goggles (1 pair)

Paper towels or other absorbent pads

Disinfectant of your choice (labeled with active ingredient and concentration e.g. 10% Bleach)

Red biohazard bags

Duct tape to seal the biohazard bags

Appropriate Personal Protective Equipment (PPE)